

ASLP™-3.5 Calcaneal Locking Plate



Indication:

- Intended use for complex fracture of calcaneus.
- Osteotomies of the calcaneus.
- For comminuted fractures.

Features:

- Anatomically designed plate matches the lateral aspect of calcaneus.
- Low implant profile reduces soft tissue irritation.
- Screw holes can accommodate 3.5 mm locking, 3.5 mm cortex and 4.0 mm cancellous screws.

ASLP-3.5 Calcaneal locking plate

PT137.60.01
| Right, L-60 mm

PT137.60.02
| Left, L-60 mm

PT137.68.01
| Right, L-68 mm

PT137.68.02
| Left, L-68 mm