

# MS - III MIS System



## Features -

- Can be performed **one-two** level of Lumbar fixation
- Preserving the posterior spinal structure with minimal invasive procedure; a small incision and marginal soft tissue dissection.
- Break-off inner nut ensures the adequate final locking of the pedicle screw
- Unique and user-friendly instrumentation ensures smooth OR procedures

# MS - III

## MIS System



### Size

MIS Poly Axial Pedicle screw	Dia 5.5 mm	Dia 6.5 mm
Length 35 mm	ST123.55.35	ST123.65.35
Length 40 mm	ST123 .55.40	ST123.6 5.40
Length 45 mm	ST123 .55.4 5	ST123.6 5.45
Length 50 mm	ST123 .55.50	ST123.6 5.50

5.5 mm Pre-contoured Rod in mm	
RT131.55.30	30
RT131.55.40	40
RT131.55.50	50
RT131.55.60	60
RT131.55.70	70
RT131.55.80	80
RT131.55.90	90
RT131.55.100	100
RT131.55.110	110
RT131.55.120	120
RT131.55.130	130
RT131.55.140	140
RT131.55.150	150

MIS INNER	
IT133	MIS Inner