



PFAN



- Proximal Diameter - 15 mm
- Medial - Lateral Angle - 4°
- CCD - Angle - 130°

Proximal Femur Antirotation Nailing System

Indication

PFAN Short :Petrochantric, Intertrochantric and High subtrochantric fractures.

PFAN Long :Low and extended subtrochantric, Combination & Pathological fractures.

PFAN Blade



Higher level of rotational & angular stability achieved through PFAN blade.

Insertion of PFAN Blade



Locked PFAN Blade.



Unlocked PFAN Blade.

Proximal Locking Jig



Distal Locking Jig




Blade Nail




Locking Screw




End Cap



End Cap II






Long PFAN



Short PFAN



Instrument Set



Miraclus Orthotech Pvt. Ltd.

Corporate office :

B-513, Mondeal Square, Opp. Karnavati Club,
Pralahadnagar, Ahmedabad-380015
Phone : +91 79 4004 8577

Plant address :

Plot No. 1112, Phase-3, G.I.D.C., Vatva, Ahmedabad - 382 445.
Phone : +91 79 2583 1332, 2583 1333
Email : contact@miraclus.com | www.miraclus.com